Manhattan Classic Gift Leotards

*Enter all leotard sizes for your gymnasts on your Meet Maker or All Gymnastics Roster using the Destira sizing guide below!

Girls' Leotards

Chest	Waist	Hip	Torso	Comfortable Fit
22 - 23"	17 1/2 - 19"	20 1/2 - 21"	38 - 40"	Child XXS (3-4)
23 - 25"	18 1/2 - 20"	21 1/2 - 23"	40 - 43"	Child XS (5-6)
25 - 26"	20 - 21 1/2"	23 - 24 1/2"	43 1/2 - 45"	Child S (6x-7)
26 - 27 1/2"	21 - 22 1/2"	24 - 25 1/2"	45 1/2 - 47*	Child M (8-10)
27 - 28 1/2*	22 - 23"	25 - 26 1/2"	47 - 49"	Child L (10-12)

Women's Leotards

Chest	Waist	Hip	Torso	Comfortable Fit
28 - 30"	19 1/2 - 21"	27 1/2 - 29"	51 - 52 1/2"	Junior (12-14)
29 1/2 - 31"	21 1/2 - 23"	29 - 30 1/2"	52 1/2 - 54"	Adult XS
30 1/2" - 33"	22 1/2 - 24"	30 1/2 - 32"	54 - 55 1/2"	Adult S
32 1/2 - 33 1/2"	24 - 25 1/2"	32 - 33"	56 - 57 1/2"	Adult M
34 - 35"	25 1/2 - 26 1/2"	33 1/2 - 35"	57 - 59"	Adult L
35 - 36"	26 1/2 - 28"	35 - 36"	59 - 61"	Adult XL

How to Measure

- Chest Measurement: taken with a deep breath, arms held out at shoulder height. Go directly over shoulder blades on back.
- 2. Waist Measurement: taken at natural waistline, on-or-around belly-button
- 3. Hip Measurement: taken around the fullest part of the buttocks
- Torso Measurement: taken from the center of the shoulder, down the front, between the legs, and up to the back to the starting point.

Rules of Thumb

- Destira leotards run fairly true to size. Most often, a gymnast will wear the same size they wear in street clothes. For growing room, buy one size up.
- Metallic fabrics have less stretch and will have a tighter feel when worn. Keep fabric type in mind when choosing size.